

Transport Stretcher - P8005

A cost-effective option across all care environments



Features and benefits

Heavy-gauge, steel-reinforced frame ensures durability, even with frequent usage

The **ability to lower the stretcher** to a height of 58.4 cm results in safer patient access

Fully integrated siderails mean safer patient transport

Optional **Steering Plus™** directional system incorporates fully-locking five castor braking, providing improved steering, stopping, cornering and control

Optional **integrated IV pole transport device** and **ergonomic push handles** allow for greater manoeuvrability

Transport stretcher - P8005



Optional: push handles



Optional: integrated IV pole



Optional: Steering Plus System

Technical Specifications (standard configuration)

Safe Working Load	317.5 kg
Weight (without surface/accessories)	120.2 kg
Height adjustment	58.4 to 87.0 cm
Trendelenburg/Reverse Trendelenburg	- 18 to + 18 °
Protection class	IPX4
Deck width	
standard	66 cm
optional	76.2 cm
Width	
siderails stored	77.2 or 87.3 cm
siderails up	81.3 or 91.4 cm
Deck length	190.5 cm
Length	210.8 cm
Surface height	
standard	7.6 cm
optional	10.2 or 12.7 cm

Siderails

height	35.6 cm
length	119.4 cm

Clearance

2.9 to 8.9 cm

Castors

4 x 20.3 cm
single band

Additional options

Push handles
Steering Plus System
4 brake & steer pedals
Integrated O² tank holder / utility tray
Radiolucent head section
Upright cassette holder
Integrated IV pole - 2 or 3 stage

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN145301-02, 22 March 2016

Enhancing outcomes for patients and their caregivers:

Hill-Rom